

# ***Pichia Yeast Transformation protocol***

## **Restriction Digest**

1. Digest ~5 – 10  $\mu$  g of plasmid DNA with one of the restriction enzymes Sac I or Sal I .
2. Check a small aliquot of your digest by agarose gel electrophoresis for complete linearization.
3. If the vector is completely linearized, heat inactivate or add EDTA to stop the reaction, phenol/chloroform extract once, and ethanol precipitate using 1/10 volume 3 M sodium acetate and 2.5 volumes of 100% ethanol. Centrifuge the solution to pellet the DNA, wash the pellet with 80% ethanol, air-dry, and resuspend in 10  $\mu$  l sterile, deionized water. Use immediately or store at  $-20^{\circ}$  C.

## **Preparing *Pichia* for Electroporation**

Follow the procedure below to prepare your *Pichia pastoris* strain for electroporation.

1. Grow 5 ml of your *Pichia pastoris* strain in YPD in a 50 ml conical tube at  $30^{\circ}$  C overnight.
2. Inoculate 500 ml of fresh medium in a 2 liter flask with 0.1 – 0.5 ml of the overnight culture. Grow overnight again to an OD600 = 1.3 – 1.5.
3. Centrifuge the cells at  $1500 \times g$  for 5 minutes at  $4^{\circ}$  C. Resuspend the pellet with 500 ml of ice-cold ( $0^{\circ}$  C), sterile water.
4. Centrifuge the cells as in Step 3, then resuspend the pellet with 250 ml of ice-cold ( $0^{\circ}$  C), sterile water.
5. Centrifuge the cells as in Step 3, then resuspend the pellet in 20 ml of ice-cold ( $0^{\circ}$  C) 1 M sorbitol.
6. Centrifuge the cells as in Step 3, then resuspend the pellet in 1 ml of ice-cold 1 M sorbitol for a final volume of approximately 1.5 ml. Keep the cells on ice and use that day. Do not store cells.

## **Transformation by Electroporation**

1. Mix 80  $\mu$  l of the cells from Step 6 (above) with 5 – 10  $\mu$  g of linearized pPICZ $\alpha$  DNA (in 5 – 10  $\mu$  l sterile water) and transfer them to an ice-cold ( $0^{\circ}$  C) 0.2 cm electroporation cuvette.
2. Incubate the cuvette with the cells on ice for 5 minutes.
3. Pulse the cells according to the parameters for yeast (*Saccharomyces cerevisiae*) as suggested by the manufacturer of the specific electroporation device being used.
4. Immediately add 1 ml of ice-cold 1 M sorbitol to the cuvette. Transfer the cuvette contents to a sterile 15 ml tube.
5. Let the tube incubate at  $30^{\circ}$  C without shaking for 1 to 2 hours.
6. Spread 50 – 200  $\mu$  l each on separate, labeled YPDS plates containing the appropriate concentration of corresponding antibiotic.
7. Incubate plates for 2 to 3 days at  $30^{\circ}$  C until colonies form.
8. Pick 10 – 20 colonies and purify (streak for single colonies) on fresh YPD or YPDS plates containing the appropriate concentration of corresponding antibiotic.