

# Macaloïd/Roche method

Protocols

Februari 2008 Version **5.2**

## 1. Introduction

- This is a combination of the macaloïd method and the Roche RNA isolation kit, for isolation of total RNA for microarray analysis.

## 2. Macaloïd part of the protocol.

1. Grow cells till wanted OD600.
2. Collect 10-20 ml culture, 1 min at 10.000 rpm (RT).

Use 10-15 OD equivalents : OD 0.1 → 100ml  
OD 1.0 → 10ml

3. Discard as much medium as possible and immediately freeze in liquid nitrogen.

If there is enough time :resuspend pellet quickly in 400 µl TE (DEPC), pipet it into a 2ml screw cup and freeze tube in liquid nitrogen.

4. Store pellet at -80°C.
5. Thaw cell pellet on ice.
6. Resuspend pellet in 400 µl TE (DEPC) if not done yet.
7. Add to screw-cap tubes containing (should be prepared in advance):
  - 0.5 g glass beads
  - 50 µl 10% SDS
  - 500 µl phenol / chloroform:IAA

Premix 300 µl chloroform:IAA (24:1) and 300 µl phenol (acid phenol, RNA grade, Sigma P4682) and take 500 µl of the organic phase.

Macaloïd layer (150-175 µl – not exact as it is very viscous)

8. Place in **bead beater**. 2x pulses ('homogenize'), 1 min interval on ice.
  - *L.lactis* 2x 45"
  - *S.pneu* and *B.sub* 2x 60"
9. Centrifuge **10** min at 10.000 rpm (4°C).
10. Transfer upper phase to fresh tube, add 500 µl chloroform:IAA (24:1)
11. Centrifuge **5** min at 10.000 rpm (4°C).

---

### 3. Roche part of the protocol.

12. Transfer 500  $\mu$ l of upper phase to fresh tube, add 2 volumes (1 ml) of lysis/binding buffer (● green cap) and mix by pipetting up and down.
13. Combine filter and collection tubes and pipette 800  $\mu$ l to upper reservoir.
14. Centrifuge 15 sec at 8600 rpm.
15. Discard flow through and pipette remainder of sample (~700  $\mu$ l) to upper reservoir.
16. Centrifuge 15 sec at 8600 rpm and discard flow through.
17. Add 100  $\mu$ l DNaseI mix. (90  $\mu$ l DNase buffer and 10 $\mu$ l DNaseI).
18. Incubate 20-30 min at 15-25 °C.
19. Add 500  $\mu$ l wash buffer I (● black cap) to upper reservoir, centrifuge as before and discard flow through.
20. Add 500  $\mu$ l wash buffer II (● blue cap) and centrifuge as before.
21. Add 200  $\mu$ l wash buffer II (● blue cap) and centrifuge 2 min at max speed.
22. Discard collection tube and insert filter tube in a sterile 1.5 ml eppendorf tube.
23. Add 50  $\mu$ l elution buffer (o colorless cap) directly on filter in upper reservoir and incubate 2-10 min at RT and centrifuge 1 min at 8600 rpm.
24. Save aliquots of the RNA at -80°C.

---

### 4. Measuring the RNA concentration

Measure concentration using the Nanodrop

1 A260 unit RNA = 40  $\mu$ g/ml

**Recommended amount for labeling is 10-20  $\mu$ g**

---

### 5. Control gel (Optional)

- Rinse casting tray and the comb with 3% H<sub>2</sub>O<sub>2</sub>
  - Pipette 150  $\mu$ l 1 M Guanidine thiocyanate in a casting tray. (Fw=118,2 118 mg/ml prepare fresh).  
**This will denature the RNA and gives better separation.**
  - Add 30 ml Agarose (1.5 % in TBE) with 2  $\mu$ l EthBr (5mg/ml)/ 100 ml to the tray and mix using the comb.  
**Use the same Agarose as for DNA gels.**
  - Load the gel :
    - 1.0  $\mu$ l totRNA
    - 2.0  $\mu$ l MQ water
    - 4.0  $\mu$ l formamide
    - 1.0  $\mu$ l loading dye (for RNA)
  - Use normal DNA reference marker.
  - Run at 8 V/cm (30mA) in 1x TBE buffer (same as DNA electrophoresis buffer)
-

---

**Sometimes it is advisable to concentrate your RNA:**

add: 1/10 volume of 1M NaAc (pH 4.5) and 2.5 volume cold EtOH

large volumes: o.n. – 20 °C

small volumes : 10' – 80 °C

---

## 6. Agilent

- Dilute 1 ul of sample with 50 ul of MiliQ DEPC to get a concentration of 20-200 ng/μl
  - Use 1 ul diluted RNA sample to check quality on the Agilent BioAnalyser
  - Ratio 23S:16S:~2.0
- 

## 7. Materials.

**acid phenol**, RNA grade  
Sigma P4682

**macaloid:**

- **Suspend 2 g. macaloid in 100ml TE**
- **Boil 5' and cool to RT**
- **Sonicate, by burst, until macaloid gels**
- **Centrifuge and resuspend in 50 ml TE pH8**
- **Store at 4°C**

Bentone MA (rheological additive)  
Rheox Inc.  
Nettlehill Road, Livingston  
West Lothian  
EH54 5 DL Scotland  
UK

**Solution :** All solution should be treated with diethyl pyro-carbonate (DEPC).  
Add 100μl DEPC/100 ml solution, incubate over night at 37°C and  
autoclave for 15 minutes.

---

**Always wear gloves when handling RNA.**

---