




## Procedure:

### Critical Steps:

- Swirl agarose to make sure that all is melted (but not too vigorously to avoid bubbles)

### The following is for a 1% Gel



0.7% agarose

- Measure 0.7 g agarose and add to flask**
- Add 70 mL 1X TAE to the flask**
- Microwave the agarose for 1:00, or until all agarose is dissolved and the solution is completely clear and bubbling**
  - If the agarose is still suspended/undissolved, microwave for additional 10-20 seconds
  - CAREFUL, bottle will be extremely hot, watch out for steam coming from the bottle, hot gloves should be by the microwave
- Cover flask with saran wrap (should be on top of the microwave)**
- Obtain gel tray, box, and comb, then rinse with dH2O**
  - Comb size depends on experiment
- Place gel tray in gel box sideways – the tray’s rubber gasket should make a watertight seal against the gel box**
  - Make sure the gasket hasn’t derailed
- Insert gel comb**
- When agarose is warm to the touch, add 7  $\mu$ L SYBR green while gently swirling**
- Once SYBR green has completely mixed, pour agarose into tray slowly**
- Use the comb or a sterile pipette tip to sweep bubbles away from the important parts of the gel, then cover the gel while it cools to prevent debris from falling in**
- Add 5X loading dye to DNA samples**
  - Bear in mind the maximum volume each type of comb can hold

- When agarose has completely solidified, gently pull tray out of box, rotate tray, then re-insert tray so the wells are on the black side**
- Pour 1X TAE into the box until it just covers the lanes**
  - Pour directly onto the lanes to expunge bubbles
- Write down the order of your DNA samples, details about your DNA samples, and where the ladder(s) will go**
  - **AVOID putting ladder in the center because you want to differentiate sides**