

## **Plasmid Extraction from *E. coli* - Alkaline Lysis Method**

(by Susanne, 2010)

- Harvest 2-4 ml of cells in eppendorf (13,000rpm, 1 min) Decant supernatant (aspirate)
- Resuspend cells in 300  $\mu$ l P1 buffer to a homogenous suspension
- Add 300  $\mu$ l of lysis buffer (P2 buffer), invert about 6 times (not more!)
- Add 300  $\mu$ l K-Ac/5% formic acid (P3 buffer) and invert tube approx 6 times. Should see a precipitate form
- Spin at 13,000 rpm for 10 min then transfer supernatant into new eppendorf
- Precipitate plasmid DNA in 0.7 vol (i.e. 630  $\mu$ l) of room temperature isopropanol and invert about 6 times
- Spin at 13,000 rpm for 15mins and decant supernatant.
- Wash pellet in 70% ethanol (ca. 700  $\mu$ l) and remove supernatant, spin again if pellet becomes dislodged.
- Quick spin to remove final trace ethanol and allow pellet to air dry (approx 10-15 mins)
- Dissolve DNA in 50-100  $\mu$ l of MQ H<sub>2</sub>O (pH5.5) or 10 mM Tris/HCl (pH8.0).

## **Recipes:**

### **P1 Buffer (Recipe from Qiagen kit)** (store in fridge)

50mM Tris/HCl [pH 8]

10mM EDTA [pH 8]

Make up part of the final volume with the Tris/HCl and EDTA solutions with water.

100µg/ml DNase-free RNase (from 10 mg/ml stock)

### **Lysis Buffer (P2)** (store at RT, but only make about 10 or 20 ml as it doesn't keep forever)

0.2M NaOH

1% SDS

### **K Acetate/5% formic acid (P3)**(store at RT)

88.3g K-acetate

15ml Formic Acid

300ml volume with dH<sub>2</sub>O