

# Sequencing Preparation

<b>Rationale:</b>	
<b>Special Observations:</b>	
<b>Results:</b>	
<b>Interpretation:</b>	

**Experiment Date:**  
**Experiment Time:**

Source: [Cornell Sequencing Handbook](#)

**Primary Experimenter (contact):**  
**Other Experimenters:**

Assembled: 6/27/2012

Reagent	Details	Quantity
ddH2O (nuclease-free)		Up to 18 $\mu$ L
Primer ( <b>8 <math>\mu</math>M stock</b> )	Conc. labeled on microfuge tubes in primer box	1 $\mu$ L
DNA template		*Var.
Sequencing tube		

\***Plasmid**,  $\sim 1 \mu$ g; **PCR product**, #base pairs/5.0 = amount of PCR product in ng that we need  
 Example: 250bp PCR product.  $250\text{bp} \div 5.0 = 50\text{ng}$  of DNA + 8 pmole primer in 18ul (Note: maximum PCR product concentration is 100ng/ul)

Order #	Template	Template conc. (ng/ $\mu$ L)	Template volume ( $\mu$ L)	Primer ID
-	Plasmid ID , Sample name; Ex) p8, Gibson 1			Ex) 0033
1				
2				
3				
4				
5				
6				
7				

<b>8</b>				
<b>9</b>				
<b>10</b>				
<b>11</b>				
<b>12</b>				
<b>13</b>				
<b>14</b>				
<b>15</b>				

**NOTE: Please write calculations on back of sheet.**

### Procedure:

#### Critical Steps:

- Clearly label all sequencing tubes before beginning
- NEVER add primers directly from blue cap tubes, only the diluted 8 μM stocks

NOTE:

For primers that aren't strictly sequencing primers, stocks are 10μM. These can still be used for sequencing.

#### Label microcentrifuge tubes

- Write the order # on the top of the tube
- Write the Order name on the side of the tube
- Ex) "Gibson 1" on side, "1" on top

#### Calculate μL of DNA template to add to reach desired amount

$$\frac{\text{Desired ng DNA}}{X \text{ DNA template } \left(\frac{\text{ng}}{\mu\text{L}}\right)} = (\text{Y } \mu\text{L to add for desired ng DNA})$$

- Aim for ~1 μg for plasmid DNA
- Please write calculations on backs of the front pages of this protocol

#### Calculate μL of ddH2O to add

- (17 μL) – (Y μL to add for desired μg DNA) = (μL H2O to add)
- Please write calculations on backs of the front pages of this protocol

#### Add ddH2O, primer, and finally DNA template to tube

- Add primer from the 8 μM stock