

Experiments:

Experiment 1

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
0	40	40F	55.76	40R	58.42	909
2	42	42F	55.62	42R	59.83	861
5	45	45F	55.53	45R	58.33	846
7	47	47F	55.29	47R	60.97	930
8	48	48F	55.02	48R	58.6	966

Experiment 2

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
58	58	45F	55.29	ExoR	57.1	1545
59	59	47F	55.62	ExoR	57.1	1611

Experiment 3

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
3	43	43F	56.43	43R	58.73	783
4	44	44F	57.21	44R	58.91	969
9	49	49F	57.24	49R	61.07	855
10	50	50F	56.95	50R	58.8	831

Experiment 4

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
57	57	44F	57.21	ExoR	57.1	1650
60	60	49F	57.24	ExoR	57.1	1536
61	61	50F	56.95	ExoR	57.1	1512

Experiment 5

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
1	41	41F	58.35	41R	59.97	1107
6	46	46F	57.67	46R	61.33	966

Protocols from Kappa:

KAPA HiFi reactions must be **set up on ice**, according to the following recommendations:

Component	Final concentration	Volume in a 25 µl reaction ¹
PCR grade water	–	Up to 25.0 µl
5X KAPA HiFi Fidelity or GC Buffer ² (contains 2.0 mM Mg ²⁺ at 1X)	1X	5.0 µl
KAPA dNTP Mix (10 mM each dNTP) ³	0.3 mM each dNTP	0.75 µl
Forward primer (10 µM)	0.3 µM	0.75 µl
Reverse primer (10 µM)	0.3 µM	0.75 µl
DMSO (100%) (for amplicons with a GC content >70%)	5%	1.25 µl
Template DNA ⁴	As needed	10 – 100 ng for genomic DNA 1 – 10 ng for less complex DNA
KAPA HiFi DNA Polymerase (1 U/µl)	0.5 U/25 µl rxn	0.50 µl

¹ For smaller reaction volumes, scale all volumes down proportionately. Reaction volumes >25 µl are not recommended.

² Use the Fidelity Buffer for most reactions. The GC Buffer is specifically formulated for amplicons/templates with a high GC content and/or stable secondary structure, and is recommended when the Fidelity Buffer produces a low yield. The additives in the GC Buffer results in a two-fold decrease in fidelity (as compared to the Fidelity Buffer).

³ The performance of all high-fidelity DNA Polymerases, including KAPA HiFi, is highly dependent on the quality of dNTPs used; the presence of even small amounts of dUTP has a dramatic impact. Only use the highest quality KAPA dNTP Mix supplied in your kit.

⁴ Use 10 ng genomic DNA or 1 ng less complex DNA per 25 µl reaction as a starting point.

Step	Temperature	Time	Number of cycles
Initial denaturation ¹	95 °C	2 – 5 min	1
Denaturation¹	98 °C	20 sec	15 - 35 ⁴
Primer annealing ²	60 – 75 °C	15 sec	
Extension ³	72 °C	15 – 60 sec/kb	
Final extension	72 °C	1 – 5 min	1
Cooling	4 – 10 °C	HOLD	1

¹ KAPA HiFi Buffers have a higher salt concentration than conventional PCR buffers. Since buffer composition affects DNA melting, it is important to denature complex templates and targets with a high GC content sufficiently. **The above denaturation parameters must therefore be strictly adhered to.**

² Due to the high salt concentration of KAPA HiFi Buffers, the optimal annealing temperature for a specific primer pair is likely to be different than when used in a conventional PCR buffer. An annealing temperature of 65 °C is recommended as a starting point. Two-step cycling protocols, with a combined annealing/extension temperature in the range of 68 – 75 °C and a combined annealing/extension time of 30 sec/kb may also be used.

³ Use 15 sec per cycle for amplicons ≤1 kb and 30 – 60 sec/kb per cycle for long amplicons or to improve yields.

⁴ ≤25 cycles are recommended for most high-fidelity applications. In cases where very low template concentrations or low reaction efficiency results in low yields, 30 or 35 cycles may be performed to produce sufficient product for downstream analysis or cloning.

Experiment #1:

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
0	40	40F	55.76	40R	58.42	909
2	42	42F	55.62	42R	59.83	861
5	45	45F	55.53	45R	58.33	846
7	47	47F	55.29	47R	60.97	930
8	48	48F	55.02	48R	58.6	966

1. Dissolve colonies 40, 42, 45, 47, 48 in 25 ul nuclease-free water separately (1 colony sample in 25ul) in 1.5ml eppendorf tubes. Heat the samples for 5 min at 95C. Centrifuge at 13,000 RPM for 3 min.
2. Make master mix for: 6 x 25ul Reactions:
Combine the following in a 500ul tube: (22ul x 6 = 132ul)
97.5 ul nuclease-free H₂O
30ul HiFi Buffer (5x)
4.5ul dNTP Mix
3. Aliquot 22 ul of master mix from previous step into 5 labeled PCR tubes (0, 2, 5, 7, 8)
4. Add .75ul 10uM Primers A and B each and 1ul template into corresponding tubes:

Tube ID	Template (SIMD#)	Primer A	Primer B
0	40	40F	40R
2	42	42F	42R
5	45	45F	45R
7	47	47F	47R
8	48	48F	48R

5. Add .5ul HiFi DNA Polymerase to each tube
6. Thermocycle as follows:

Initial denaturation: 4min 30sec at 95C

35 cycles of:

Denaturation: 20sec at 98C

Annealing: 15sec at 55

Extension: 30sec at 72C

Final Extension: 2min at 72C

Experiment #2:

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
58	58	45F	55.29	ExoR	57.1	1545
59	59	47F	55.62	ExoR	57.1	1611

1. Follow same directions as in Experiment 1 for colonies 58 and 59 but make 66ul of master mix:
Combine the following in a 500ul tube: (2 x 22ul)
48.75ul nuclease-free H₂O
15ul HiFi Buffer (5x)
2.25ul dNTP Mix

3. Aliquot 22 ul of master mix from previous step into 2 labeled PCR tubes (58, 59)
4. Add .75ul 10uM Primers A and B each and 1ul template into corresponding tubes:
5. Add .5ul HiFi DNA Polymerase to each tube
6. Thermocycle as follows:

Initial denaturation: 4min 30sec at 95C

35 cycles of:

Denaturation: 20sec at 98C

Annealing: 15sec at 55

Extension: 1min 20sec at 72C

Final Extension: 3min at 72C

Experiment #3:

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
3	43	43F	56.43	43R	58.73	783
4	44	44F	57.21	44R	58.91	969
9	49	49F	57.24	49R	61.07	855
10	50	50F	56.95	50R	58.8	831

2. Dissolve colonies 43, 44, 49, and 50 in 25 ul nuclease-free water separately (1 colony sample in 25ul) in 1.5ml eppendorf tubes. Heat the samples for 5 min at 95C. Centrifuge at 13,000 RPM for 3 min.

2. Make master mix for: 5 x 25ul Reactions:

Combine the following in a 500ul tube: (22ul x 5 = 110ul)

81.25ul nuclease-free H2O

25ul HiFi Buffer (5x)

3.75ul dNTP Mix

3. Aliquot 22 ul of master mix from previous step into 5 labeled PCR tubes (3,4,9,10)
4. Add .75ul 10uM Primers A and B each and 1ul template into corresponding tubes:
5. Add .5ul HiFi DNA Polymerase to each tube
6. Thermocycle as follows:

Initial denaturation: 4min 30sec at 95C

35 cycles of:

Denaturation: 20sec at 98C

Annealing: 15sec at 56C

Extension: 30sec at 72C

Final Extension: 2min at 72C

Experiment #4:

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
57	57	44F	57.21	ExoR	57.1	1650
60	60	49F	57.24	ExoR	57.1	1536
61	61	50F	56.95	ExoR	57.1	1512

Experiment #5:

Rxn ID	Strain: SIMD#	Primer A	Primer A TM	Primer B	Primer B TM	Expected Product Size in bp
1	41	41F	58.35	41R	59.97	1107
6	46	46F	57.67	46R	61.33	966

3. Dissolve colonies 41, 46, 57, 60, 61 in 25 ul nuclease-free water separately (1 colony sample in 25ul) in 1.5ml eppendorf tubes. Heat the samples for 5 min at 95C. Centrifuge at 13,000 RPM for 3 min.

2. Make master mix for: 6 x 25ul Reactions:

Combine the following in a 500ul tube: (22ul x 6 = 132ul)

97.5 ul nuclease-free H₂O

30ul HiFi Buffer (5x)

4.5ul dNTP Mix

3. Aliquot 22 ul of master mix from previous step into 5 labeled PCR tubes (1, 6, 57, 60, 61)

4. Add .75ul 10uM Primers A and B each and 1ul template into corresponding tubes:

5. Add .5ul HiFi DNA Polymerase to each tube

6. Thermocycle as follows:

Experiment 4 samples:

Initial denaturation: 4min 30sec at 95C

35 cycles of:

Denaturation: 20sec at 98C

Annealing: 15sec at 56.5C

Extension: 1min 20sec at 72C

Final Extension: 3min at 72C

Experiment 5 samples:

Initial denaturation: 4min 30sec at 95C

35 cycles of:

Denaturation: 20sec at 98C

Annealing: 15sec at 57C

Extension: 1min at 72C

Final Extension: 2min 30sec at 72C