

Vent PCR

Rationale:	
Special Observations:	
Results:	
Interpretation:	

Experiment Date:

Source: [NEB](#)

Experiment Time:

Primary Experimenter (contact):

Assembled: 6/27/2012

Other Experimenters:

Reagent	Details	Quantity	
		Suggested:	Used:
ddH ₂ O (nuclease-free)		*Var.	
dNTP mix (10 mM)		1 µL	
10X ThermoPol Reaction Buffer		5 µL	
Forward Primer (10 µM)	(ID)→	1 µL	
Reverse Primer (10 µM)	(ID)→	1 µL	
Template DNA	(Name)→	**1 µL	
DMSO (optional)	Prevents primer secondary structure, alter H ₂ O if	3%	
Vent DNA polymerase		0.5 µL	
MgSO ₄		***1-6mM	
		50 µL Total	

*H₂O volume to use for non-DMSO rxn = (50 µL) – (9.5 µL) – (µL MgSO₄)

**1 µl of 1-10 ng/µl for vector DNA; 1 µl of 50-250 ng/µL for genomic DNA

*** Primer extensions longer than 2 kb almost always require magnesium levels higher than 2 mM, while for primer extensions shorter than 2 kb, there is no correlation between length and optimum magnesium concentration

Critical Steps:

- Add Phusion last, minimize time out of freezer
- Program PCR machine before adding Phusion
- Add each component in order listed above to a PCR tube

NOTE:

- Make a mastermix for number of PCRs + 1 if doing more than two PCRs (mastermix includes H₂O, dNTPs, PCR buffer)
- See: [NEB Vent PCR optimization](#) if issues arise

PCR Machine Settings:

		Recommended:	Used:	Recommended:	Used:
Step 1	Initial denaturing	95 °C		5 minutes	
Step 2 (25 – 30 cycles)	Denature	95 °C		30 seconds	
	Anneal	*Variable		30 seconds	
	Extend	72 °C		60 sec/kb	
Step 3	Final Extension	72 °C		5 minutes	
Step 4	Hold	4 °C		Indefinite	

* Lowest primer Tm – 5 °C