

Quantitative real-time PCR (qPCR)

Sample preparation:

- 5 µL primary virus stock solution
- 7,5 µL DNaseI
- 5 µL MgCl₂ [50 mM]
- add 50 µL RNase free water

→ 30 min incubation at 37 °C

→ 10 min enzyme inactivation at 65 °C

Standard curve preparation:

1. Calculation of the plasmid weight

$$m_{\text{plasmid}} = n \text{ [bp]} * 1 \text{ mole} / 6.096 * 10^{23} * 660 \text{ g/mole}$$

$$m_{\text{plasmid}} = n \text{ [bp]} * 1.096 * 10^{-21} \text{ [g/bp]}$$

2. Calculation of the copies of the plasmid

$$m_{\text{copies}} = m_{\text{plasmid}} * \text{copies}$$

→ calculate copies from $1.3 * 10^1$ to $1.3 * 10^7$ copies

QPCR preparation:

Pipetting scheme for qPCR:

component	Volume per reaction [µL]
2 x QuantiFast SYBR Green pcr mix	7.5
Forward primer	1 [10 µM]
Reverse primer	1 [10 µM]
Template	5.0
Rnase free water	0.5
Total volumen	15

- preparing duplicates of standard curve sample and negative control

QPCR protocol:

- two step cycling program

Step	Time	Temperature [°C]
Initial activation	5 min	95
Denaturation	10 s	95
Annealing/extension	30 s	60

Number of cycles: 45