

# Site Directed Mutagenesis (QuickChange Method)

## Cornell iGEM 2012 Protocol

Source: Dylan Webster (Adapted from QuickChange II XL Site-Directed Mutagenesis Kit Protocol)

<http://www.chem.agilent.com/Library/usermanuals/Public/200521.pdf>

<http://www.neb.com/nebecomm/products/protocolProductE0553.asp>

### Mutagenic Primer Design:

Use the Agilent webpage to help you design mutagenic primers: [www.agilent.com/genomics/qcpd](http://www.agilent.com/genomics/qcpd)

In general, primer pairs should anneal to the same sequence of opposite strands, contain the desired mutation flanked by 10-15 bases of correct sequence on both sides, and have a melting temperature  $\geq 78^\circ\text{C}$ .

### Mutant Strand Synthesis Reaction:

\* Ensure that miniprep plasmid came from a dam<sup>+</sup> strain (DH5 $\alpha$  is good).

Prepare the following 50  $\mu\text{L}$  reaction **on ice**:

- 10  $\mu\text{L}$  of 5X Phusion HF or GC Buffer (Try HF first).
- 1  $\mu\text{L}$  10 mM dNTPs.
- 2.5  $\mu\text{L}$  sense primer.
- 2.5  $\mu\text{L}$  antisense primer.
- 1.5  $\mu\text{L}$  DMSO.
- X  $\mu\text{L}$  template DNA (Try 10ng first).
- X  $\mu\text{L}$  ddH<sub>2</sub>O (Up to 50  $\mu\text{L}$ )
- 0.5  $\mu\text{L}$  Phusion DNA polymerase (1.0 unit, **ADD LAST**).

Gently flick the PCR tube to mix the reaction mixture, and perform a quick spin in the Bio-Rad microcentrifuge. Put the PCR tube back on ice after spinning down, and preheat the thermal cycler to  $98^\circ\text{C}$  before transferring the PCR tube. Cycle the reaction according to the following parameters:

Segment	Cycles	Temperature	Used:	Time	Used:
1	1	$98^\circ\text{C}$		1 minute	
2	18	$98^\circ\text{C}$		30 seconds	
		* $60-65^\circ\text{C}$		50 seconds	
		$72^\circ\text{C}$		15-30 sec/kb	
3	1	$72^\circ\text{C}$		10 minutes	
4	1	$4^\circ\text{C}$		Indefinite	

\* This is experimental, since the polymerase that the QuickChange kit uses tends to need a lower annealing temperature than our Phusion does. The QuickChange protocol calls for a  $60^\circ\text{C}$  annealing temperature. We may need to optimize this empirically.

**Make sure the PCR reaction has been at  $4^\circ\text{C}$  for at least 5 minutes before proceeding. Also, make sure to have a  $37^\circ\text{C}$  incubator ready to go!**

*Optional:* 10  $\mu\text{L}$  of PCR reaction can be checked on a gel. QuickChange protocol says to not worry about it if bands aren't visualized at this stage, and to proceed with DpnI digestion in either case.

### DpnI Digestion of the PCR Products

- Add 1  $\mu$ L DpnI directly to the PCR reaction in an ice bucket.
- Pipette up and down and flick to mix the reaction, and then spin down the mixture in the Bio-Rad microcentrifuge.
- As quickly as possible, transfer the reaction to a 37°C incubator (either the thermal cycler, the water bath, or the heat block).
- Heat inactivate DpnI by incubating at 80°C for 20 minutes. (The easiest way to do this is to set up a program on the thermal cycler for both incubations.) [*This step is not necessary if DNA is to be spin-column purified prior to ligation.*]

### Electroporation of Mutated Plasmid into DH5 $\alpha$

- Prior to electroporation, it is necessary to desalt the reaction mixture. This can be done either via membrane dialysis or via spin-column purification. If the spin-column method is used, I (Dylan) think that DpnI is more likely to be excluded from the desalted plasmid. I'm not sure if this will greatly enhance transformation efficiency, but it's possible.
- *Possibility to loop back up to another mutant strand synthesis reaction from here to perform multiple mutations on the same template. I'd quantify the cleaned-up DNA before doing another PCR. There'd also be enough DNA to loop back AND proceed with a transformation just in case the loop method doesn't work.*
- Thaw cells on ice while chilling electroporation cuvettes.
- Add 2  $\mu$ L desalted DNA to thawed cells, gently mix, and keep on ice for a few minutes.
- Transfer mixture to the gap of a chilled cuvette.
- Pulse at 1.8 kV, etc.
- Immediately add 950  $\mu$ L sterile SOB or SOC (heated to 37°C) to cuvette and resuspend cells.
- Transfer to a culture tube and recover in the shake incubator for 1 hour.
- Plate on antibiotic plate and incubate for about 16 hours.