

Single Digest

Rationale:	
Special Observations:	
Results:	
Interpretation:	

Experiment Date:

Source: NEB

Experiment Time:

Primary Experimenter (contact):

Assembled: 6/27/2012

Other Experimenters:

Reagent	Details	Quantity
Sterile H ₂ O		Up to 50 μ L
*10X NEB Buffer	See: Enzyme Chart to choose buffer	5 μ L
**100X BSA	See: Enzyme Chart to decide if needed	0.5 μ L
1-10 μ g DNA (Or 200 ng for minimal gel visualization)	(name)→	Var.
Restriction enzyme	(enzyme)→	1 μ L

Procedure:

Critical Steps:

- Restriction enzymes are expensive! Leave frozen until final step.
- Use small volume tubes
- Carefully label tubes
- All steps on ice
- See: Enzyme Chart to choose reaction temperature
-

NOTE:

- BSA does not inhibit any restriction enzyme

Turn on water bath

- Check enzyme chart for reaction temperature

Calculate DNA volume to use

- $(? \mu\text{L DNA}) = \frac{1000 \text{ ng}}{\text{DNA sample concentration} \frac{\text{ng}}{\mu\text{L}}}$

Calculate H2O volume to use

- $(? \mu\text{L H}_2\text{O}) = 50 - (? \mu\text{L DNA}) - 6 \mu\text{L} - (0.5 \mu\text{L if using BSA})$

Add (? $\mu\text{L H}_2\text{O}$) to reaction tube

Add 5 μL 10X NEB buffer to reaction tube

IF REQUIRED, add 0.5 μL 100X BSA to reaction tube

Add (? $\mu\text{L DNA}$) to reaction tube

Add 1 μL restriction enzyme to reaction tube

Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube. Follow with a quick ("touch") spin-down in a microcentrifuge.

- Do not vortex the reaction.

Incubate 1 hour in water bath

- Use optimal reaction temperature

Stop reaction

- If further manipulating DNA NOT required, do DNA gel electrophoresis with a loading dye that includes EDTA
- If further manipulation required, heat inactivate (See: Enzyme Chart)