

Picking clones

1. number clones of interest at the bottom of the plate
2. pick clones from plate
3. inoculate clone in 5 ml LB-medium
4. over-night incubation at 37 °C

Picking clones

- experiment date: ____ ; time: ____
- name of investigator: ____

- Clone ID_1: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_2: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_3: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_4: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_5: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_6: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_7: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_8: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_9: Strain:_____, Plasmid:_____, Clone number: _____
- Clone ID_10: Strain:_____, Plasmid:_____, Clone number: _____